

Patient Terminal | Patient head Set

NBW-B1



Features

Thin body, including handheld, wall-mounted installation on the Patient room

Call/cancel: support calling the nurse station.

Duplex intercom: support full duplex intercom with the host of the nurse station.

Intercom mode: adopt echo cancellation technology, support full-duplex intercom, and the intercom audio sampling rate is not less than 16KHz.

Various interfaces: support external waterproof button, door lamp.

Power supply mode: bus system power supply.

Precautions: During installation, the distance between the hand-held pager and the extension should be > 15cm.

Technical Specifications

Bus protocol : FSK, DTMF	Frequency response: 200Hz~16KHz
Harmonic distortion : ≤0.3%	MIC input sensitivity : 10mV
SNR : > 70dB	Work temperature : -10~+60 degrees
Line Out output impedance : 300Ω	Working humidity : 10%~90%
Line In input level : 16V~36V	Power : ≤2W
Signal line impedance: 300 Ω	
Built-in speaker impedance and power : 20 Ω , 1W	
Hand-held length: 2m (maximum tension); 1.2m (hanging)	
Size of Terminal: 78x120x12mm	
Size of hand-held: 105x35x15mm	

Dimension



Position of assembly

